



productivity

ddRFormula® Series

beyond conventional DR

Swissray is the pioneer and worldwide leader in the design, manufacturing and marketing of state-of-the-art Digital Radiography systems. Swissray's unique APS™ – Automated Positioning System – takes performance beyond conventional DR by streamlining the radiography workflow to the maximum. With an industry leading 3.5 lp/mm, Swissray delivers unmatched image quality at lowest radiation dose.



Variable SID (FFD)

LCD Panel with Patient and Positioning Information

ddArt™ Backlit Cover

C-Arm Design

FP-5000™ Si Flat Panel Detector
with 3P™ – Panel Protection Program

«The C-arm design from Swissray has been very helpful in positioning our pediatric orthopedic patients. They don't seem to be threatened at all by the equipment and it's easy to get those difficult views that are hard to convince a child, are important.»

John Strain MD, Dept. Chairman, The Children's Hospital, Denver CO, USA

Ultimate **workflow**

The **ddRFormula**[®] is Swissray's innovative C-arm designed digital radiography system. Featuring the exclusive APS[™] – Automated Positioning System – it is the most automated DR solution in the marketplace. All system movements are motorized and software controlled resulting in total automatic functionality. The **ddRFormula**[®] features a C-arm design with the X-ray tube always centered to the detector for fastest, most precise and convenient patient positioning. The **ddRFormula**[®] Plus version allows for convenient off-detector imaging.

Features

- APS[™] Automated Positioning System
- Minimum Time Between Exposure («MTBE») of 3 seconds
- Exposure to Diagnostic Image («EDI») of 5 seconds
- Performs the workload of four cassette based systems
- AutoStitching for Orthopedic Imaging
- Off-center Imaging
- EASI Function – Examination Automated Sequence Imaging



from RIS/HIS



eXpert™

APS™

SwissVision®



to PACS



APS™ – Automated Positioning System

Automatically **faster**

Swissray's unique APS™ – Automated Positioning System – streamlines the radiography workflow process by automating all positioning and image acquisition requirements. Patient data can be transferred directly from the RIS/HIS via DICOM worklist, while all exposure and image processing parameters can be chosen with simple touch-screen selections. Swissray's advanced robotics position the system for the selected examination by remote control while an integrated video camera monitors the patient to ensure correct positioning.

Features

- Handheld remote controlled system positioning
- 944 pre-programmed, fully customized positions
- Advanced robotics for selected examination
- Integrated video camera to monitor patient and ensure positioning
- Memory function for consistent follow-up examinations
- Virtually eliminates repeats and retakes





FP-5000™ Detector

More **detail**, less dose

With 3.5 lp/mm spatial resolution, highest DQE and the five field measuring chamber, Swissray's FP-5000™ Silicon (Si) Flat Panel detector delivers the best diagnostic image quality at lowest radiation dose. It combines amorphous silicon technology with a cesium iodide scintillator, providing highest detail detectability with 143 micron pixels. The unique 3P™ – Panel Protection Program – provides superb collision protection thus eliminating concerns over reliability of flat panel technology.

Features

- 3.5 lp/mm spatial resolution
- 143 micron pixels
- 17" (43 cm) square format for large field-requirements
- Highest DQE
- Carbon fiber casing lowers radiation dose requirement
- Five field measuring chamber
- Featuring 3P™ with a built-in infrared sensor to detect obstacles and a floating frame shock-absorber for maximum panel protection



Ultimate panel protection in case of a collision – the 3P™ floating frame shock-absorber and an energy dampening mechanism absorb and dissipate impact energy.



The infrared sensor detects obstacles and automatically slows down and halts the system. Swissray's grid can be conveniently attached and removed in any detector position.



Patient Data Management

Connecting seamlessly

The **eXpert™ control desk** automates every aspect of the radiographic procedure. Patient demographic data is transferred directly from RIS/HIS via DICOM worklist while all exposure and image processing parameters can be chosen with a few touch screen selections.

The **SwissVision® modality workstation** comprehensively manages all patient data. Open system architecture and «IHE» proven DICOM 3.0 compliant interfaces seamlessly integrate with existing and future network connections like PACS, RIS, local workstations and modality archives.

eXpert™ features

- MPPS (Modality Performed Procedure Step): work orders are received directly from RIS/HIS via DICOM.
- Protocol driven: individual parameter preferences can be stored for multiple users
- EASI Function – Examination Automated Sequence Imaging

SwissVision® features

- Functions such as windowing, leveling, zoom, rotation and positive/negative displaying
- TrueSize Imaging
- Automatic algorithm selection minimizes post-processing requirements
- Patient data stored in DICOM header for future examinations
- Repeat/reject examination analysis
- Exposure index to monitor image quality in relation to radiation dose
- HL7 interface capability for full integration with HIS



The eXpert™ control desk with the built-in eXpertView™ camera allows for utmost positioning accuracy.



The CutOff and SizeWise functions reduce image file size and optimize data storage capacity. Text annotations can be added for convenient image identification.



Advanced Imaging Applications

Specific **solutions**

Swissray offers a wide variety of specific solutions to meet the diverse requirements of radiography.

Pediatric Imaging Solutions. Children are among the most important beneficiaries of Swissray's lowest radiation dose. The Pediatric Package includes 3.5 lp/mm spatial resolution with special imaging algorithms and X-ray parameters providing unrivalled image quality with highest detail detectability. With the ability to perform the examination in a fraction of the time, Swissray has minimized much of the fear and anxiety associated with Pediatric imaging.


Chest Imaging. Special chest imaging parameters and algorithms make the *ddRFormula*® Series the diagnostic tool of choice for pulmonary disease detection. Healthcare systems around the globe have chosen Swissray DR systems for their TB screening programs.



The *ddRFormula*® Series perform off-center imaging for orthopedic and pediatric applications.



Special chest imaging algorithms provide unmatched diagnostic information.



AutoStitching combines up to four images.

Orthopedic Excellence

Simply more **efficient**

AutoStitching. The *ddRFormula*[®] systems perform full body imaging with their unique «AutoStitching» function, automatically combining up to four images. Orthopedic studies such as Scoliosis and Long Leg imaging are performed with greater speed and precision than ever before. The EasyCorrect tool allows the user to ensure perfect alignment even in case of patient movement.

Bone Mineral Density Analysis. Software algorithms automatically analyze a digital hand radiograph to measure bone mineral density with accuracy and precision. Providing unsurpassed access, convenience and clinical performance, the analysis costs only a fraction of other BMD testing devices available. Swissray's BMD Analysis Package has its own reimbursement code and provides additional revenues.



Scoliosis Imaging



Long Leg Imaging



Bone Mineral Density Analysis



ddArt™

Creating **ambience**

ddArt™ is an exclusive backlit design on the front cover of the system. Artwork such as sports themes, cartoons or institutional banners can be selected or customized. ddArt™ greatly enhances the ambience of the radiographic examination room, thus creating a positive patient experience.



Motorized X-ray tube rotation enables off-detector imaging on cassette based media for the occasional patient that cannot be transferred to the imaging table. This feature makes the *ddRFormula*® Plus ideal for Emergency Room applications.



ddRFormula® Applications

Providing **productivity**



beyondDR

Swissray Medical AG reserves the right to make changes in specifications and/or to discontinue any product at any time without prior notice or obligation and will not be liable for any consequences resulting from the use of this publication. Any technical data contained in this document may vary within defined tolerances. Original images always lose a certain amount of detail when reproduced.

© 2007 Swissray Medical AG
All rights are reserved. Reproduction in whole or in part is prohibited without the prior written consent of the copyright holder.
SHM_35_700_001_04_E 11/07

Manufactured by

Swissray Medical AG
Turbistrasse 25–27
CH-6280 Hochdorf
Switzerland
Phone +41 41 914 12 12
Fax +41 41 914 12 13
info@swissray.com
www.swissray.com

Distributed by

Swissray International, Inc.
One Tower Center Blvd.
East Brunswick, NJ 08816
USA
Phone 908 353 0971
Toll Free 800 903 5543
Fax 908 353 1237
info@swissray.com

Swissray 